


# Pairing Your MemorialCare Omron EVOLV<sup>®</sup> – Blood Pressure Monitor

Only pair your MemorialCare blood pressure monitor after you've successfully connected your MemorialCare Virtual Exam Kit. You will only need to perform this setup once.



**Step 1:** Activate Bluetooth<sup>®</sup> on the mobile device or tablet where the TytoCare app is installed.

**Step 2:** Log in to myChart and launch the TytoCare app by selecting TytoCare Pairing & Resources in the myChart menu.

**Step 3:** In the TytoCare app, access the user menu from the top left and select **Pair External Devices**. Click on **BP Monitor** and follow the steps to pair your device.

The Bluetooth<sup>®</sup> word mark and logos are registered trademarks owned by Bluetooth SIG, Inc.